

Anti-Cd1D Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - P15813

Item Number: CM0004761



Introduction

Recombinant rabbit monoclonal antibody against human CD1d, suitable for WB, IF/ICC, FC, ELISA. Targets antigen-presenting glycoprotein involved in glycolipid presentation to NKT cells. Validated for human samples.

[Learn More](#)

Parameter	Value
Product Name	CD1d Rabbit mAb
Remarks/Alias	R3; CD1A; R3G1; CD1d
Species	Human
GeneID (human)	912
GeneID	912
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 20-301 of human CD1d (NP_001757.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	P15813
Protein Weight	38kDa
Shipping	Ice bag