


Anti-Cd1D Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P11609 / P11610

Item Number: CM0004196



Introduction

High-quality Anti-CD1d rabbit monoclonal antibody validated for WB, IF-P, and ELISA. Reacts with mouse and rat. Suitable for immunological studies of glycolipid antigen presentation and NKT cell function.

[Learn More](#)

Parameter	Value
Product Name	CD1d Rabbit mAb
Remarks/Alias	Cd1a; Cd1d; CD1.1; Ly-38
Species	Mouse
Gene ID (Human)	912
Gene ID	12479/12480
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-305 of mouse CD1d (NP_031665.2/NP_031666.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P11609/P11610
Protein Weight	39 kDa/38 kDa
Shipping	Ice bag