


Anti-Cd1D Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P15813

Item Number: CM0019312



Introduction

Human CD1D rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Cross-reacts with mouse. A useful tool for studying glycolipid presentation and NKT cell biology (SWISS P15813).

[Learn More](#)

Parameter	Value
Product Name	CD1D Rabbit pAb
Remarks/Alias	R3; CD1A; R3G1; CD1D
Species	Human
GeneID (Human)	912
GeneID	912
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-190 of human CD1D (NP_001757.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P15813
Protein Weight	38kDa
Shipping	Ice bag