

Anti-Cd2 Rabbit Monoclonal Antibody For Wb, Fc, Elisa - P08920

Item Number: CM0007029

Introduction

CD2 rabbit monoclonal antibody validated for WB, FC, ELISA with mouse cross-reactivity. Recombinant protein immunogen targeting extracellular domain. Ideal for T-cell adhesion and signaling research.

[Learn More](#)

Parameter	Value
Product Name	CD2 Rabbit mAb
Remarks/Alias	Ly37; LFA-2; Ly-37; CD2
Species	Mouse
Gene ID (human)	914
Gene ID	12481
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 23-203 of mouse CD2 (NP_038514.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Mouse
SWISS	P08920
Protein Weight	38kDa
Shipping	Ice bag