

Anti-Cd2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P06729

Item Number: CM0006011



Introduction

CD2 Rabbit mAb (P06729) is a validated rabbit monoclonal antibody for WB, IHC-P, and ELISA, targeting human CD2, a key T-cell adhesion molecule. Ideal for immunology and diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	CD2 Rabbit mAb
Remarks/Alias	T11, SRBC, LFA-2, CD2
Species	Human
GeneID (human)	914
GeneID	914
Immunogen	Recombinant protein of human CD2
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P06729
Protein Weight	39kDa
Shipping	Ice bag