

Anti-Cd2 Rabbit Polymab® Monoclonal Antibody For Wb, Ihc-P, Elisa - P06729/P08920

Item Number: CM0004787



Introduction

CD2 Rabbit PolymAb® monoclonal antibody, validated for WB, IHC-P, ELISA. Recognizes human and mouse CD2 (T11, SRBC, LFA-2), a T-cell adhesion and signaling molecule. Ideal for immunology studies.

[Learn More](#)

Parameter	Value
Product Name	CD2 Rabbit PolymAb®
Remarks/Alias	T11; SRBC; LFA-2
Species	Human/Mouse
GeneID (Human)	914
GeneID	914/12481
Immunogen	Recombinant protein of human/mouse CD2.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P06729/P08920
Protein Weight	38kDa/39kDa
Shipping	Ice bag