


Anti-Cd20 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P11836

Item Number: CM0021090



Introduction

CD20 Rabbit pAb raised against human CD20 intracellular region; validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human and mouse. Ideal for B-lymphocyte and lymphoma research.

[Learn More](#)

Parameter	Value
Product Name	CD20 Rabbit pAb
Remarks/Alias	B1; S7; Bp35; CD20; FMC7; CVID5; LEU-16
Species	Human
Gene ID (Human)	931
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P11836
Protein Weight	33kDa
Shipping	Ice bag