

Anti-Cd20 Rabbit Monoclonal Antibody For Wb And Ihc-P - P11836

Item Number: CM0005721



Introduction

Rabbit monoclonal antibody against human CD20, suitable for WB, IHC-P, IF/ICC, FC, ELISA. Validated in human and cynomolgus monkey. Targets B-lymphocyte antigen CD20, involved in B-cell activation and calcium influx.

[Learn More](#)

Parameter	Value
Product Name	CD20 Rabbit mAb
Remarks/Alias	B1; S7; Bp35; CD20; FMC7; CVID5; LEU-16
Species	Human
GeneID (Human)	931
GeneID	931
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 100-200 of human CD20 (NP_068769.2))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, FC, ELISA
Cross Reactivity	Human, Cynomolgus monkey
SWISS	P11836
Protein Weight	33kDa
Shipping	Ice bag