


Anti-Cd200 Polyclonal Antibody For Wb, Elisa - 054901

Item Number: CM0020717



Introduction

Rabbit polyclonal antibody against mouse CD200 (OX-2 membrane glycoprotein), validated for Western blot and ELISA, cross-reacts with rat. Suitable for T-cell costimulation and myeloid cell regulation research.

[Learn More](#)

Parameter	Specification
Product Name	CD200 Rabbit pAb
Remarks/Alias	OX2; Mox2; CD200
Species	Mouse
Gene ID	17470
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-232 of mouse CD200 (NP_034948.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O54901
Protein Weight	31kDa
Shipping	Ice bag