

Anti-Cd200 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P41217

Item Number: CM0007002

Introduction

Recombinant rabbit monoclonal antibody targeting human CD200 (OX-2 membrane glycoprotein). Validated for WB, IHC-P, and ELISA; cross-reacts with mouse and rat. Ideal for immunology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	CD200 Rabbit mAb
Remarks/Alias	MRC; MOX1; MOX2; OX-2; CD200
Species	Human
GenelD (human)	4345
GenelD	4345
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-232 of human CD200. (NP_001352780.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41217
Protein Weight	31kDa
Shipping	Ice bag