


Anti-Cd200R1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Td46

Item Number: CM0021025



Introduction

Human CD200R1 rabbit polyclonal antibody validated for WB and ELISA. Inhibitory receptor limiting inflammation, expressed on immune cells. Bulk IVD raw material supply from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	CD200R1 Rabbit pAb
Remarks/Alias	OX2R; MOX2R; CD200R; HCRTR2; CD200R1
Species	Human
Gene ID (Human)	131450
Gene ID	131450
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-240 of human CD200R1 (NP_740750.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8TD46
Protein Weight	37kDa
Shipping	Ice bag