


## Anti-Cd204 Monoclonal Antibody For Flow Cytometry - P21757

Item Number: CM0005634



### Introduction

Rabbit monoclonal -conjugated antibody targeting human CD204 (MSR1), validated for flow cytometry. Detects macrophage scavenger receptor, ideal for atherosclerosis and innate immunity studies.

[Learn More](#)

| Parameter        | Value                                                          |
|------------------|----------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD204 mAb                                    |
| Remarks/Alias    | SRA; SR-A; CD204; SR-AI; phSR1; phSR2; SCARA1; SR-AII; SR-AIII |
| Species          | Human                                                          |
| Gene ID (Human)  | 4481                                                           |
| Gene ID          | 4481                                                           |
| Immunogen        | Recombinant protein                                            |
| Source           | Rabbit                                                         |
| Category         | Monoclonal Antibodies                                          |
| Application      | FC                                                             |
| Cross Reactivity | Human                                                          |
| SWISS            | P21757                                                         |
| Protein Weight   | 40kDa/43kDa/50kDa                                              |
| Shipping         | Ice bag                                                        |