


Anti-Cd204 (Msr1/Scara1) Rabbit Monoclonal Antibody For IHC-P, Fc, Elisa - P21757

Item Number: CM0007345



Introduction

Anti-CD204 Rabbit mAb validated for IHC-P, FC, ELISA. Recognizes human macrophage scavenger receptor types I and II (MSR1/SCARA1), key in LDL endocytosis and foam cell formation. Suitable for atherosclerosis and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	CD204 Rabbit mAb
Remarks/Alias	SRA; SR-A; CD204; SR-AI; phSR1; phSR2; SCARA1; SR-AII; SR-AIII
Species	Human
Gene ID	4481 (Human)
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 77-451 of human CD204 (NP_619729.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, FC, ELISA
Cross Reactivity	Human
SWISS	P21757
Protein Weight	40kDa/43kDa/50kDa
Shipping	Ice bag