


Anti-Cd205/Dec-205 Rabbit Monoclonal Antibody - O60449

Item Number: CM0005717



Introduction

Rabbit monoclonal antibody against human CD205/DEC-205, validated for WB, IF/ICC, and ELISA. The endocytic receptor involved in antigen processing; suitable for immunology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CD205/DEC-205 Rabbit mAb
Remarks/Alias	CD205; LY-75; CLEC13B; DEC-205; GP200-MR6
Species	Human
GeneID (Human)	4065
GeneID	4065
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 216-501 of human CD205/DEC-205 (NP_002340.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O60449
Protein Weight	198kDa
Shipping	Ice bag