


Anti-Cd209 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Nnx6

Item Number: CM0006142



Introduction

Rabbit monoclonal antibody against human CD209 (DC-SIGN) for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying dendritic cell function and pathogen recognition.

[Learn More](#)

Parameter	Value
Product Name	CD209 Rabbit mAb
Remarks/Alias	CDSIGN; CLEC4L; DC-SIGN; DC-SIGN1; hDC-SIGN; CD209
Species	Human
GeneID (Human)	30835
GeneID	30835
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CD209 (Q9NNX6)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NNX6
Protein Weight	46 kDa
Shipping	Ice bag