


Anti-Cd209 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nnx6

Item Number: CM0015516



Introduction

CD209 (DC-SIGN) rabbit polyclonal antibody for WB, IF/ICC, ELISA. Recognizes human, mouse, rat. Recombinant immunogen. For studies of dendritic cell-mediated immunity and pathogen attachment.

[Learn More](#)

Parameter	Value
Product Name	CD209 Rabbit pAb
Remarks/Alias	CDSIGN; CLEC4L; DC-SIGN; DC-SIGN1; hDC-SIGN; CD209
Species	Human
GeneID (Human)	30835
GeneID	30835
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 265-404 of human CD209 (NP_066978.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NNX6
Protein Weight	46kDa
Shipping	Ice bag