

Anti-CD213A2/IL-13 R Alpha 2 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q14627

Item Number: CM0009786



Introduction

Recombinant rabbit monoclonal antibody targeting human CD213a2/IL-13R alpha 2, validated for WB, IP, ELISA. Cross-reacts with rat. Ideal for studying IL-13/IL-4 signaling, immune regulation, and inflammation.

[Learn More](#)

Parameter	Value
Product Name	CD213a2/IL-13 R alpha 2 Rabbit mAb
Remarks/Alias	CT19; IL-13R; IL13BP; CD213A2
Species	Human
Gene ID (Human)	3598
Gene ID	3598
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-355 of human CD213a2/IL-13 R alpha 2 (NP_000631.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Rat
SWISS	Q14627
Protein Weight	44kDa
Shipping	Ice bag