

Anti-Cd215/IL-15R Alpha Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13261

Item Number: CM0028186



Introduction

Rabbit polyclonal antibody against human CD215/IL-15R alpha (IL-15 receptor subunit alpha). Validated for Western blot, immunofluorescence/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for studying IL-15 signaling, immune cell function, and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	CD215/IL-15R alpha Rabbit pAb
Remarks/Alias	CD215; CD215/IL-15R alpha
Species	Human
GeneID (human)	3601
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13261
Protein Weight	28kDa
Shipping	Ice bag