

Anti-Cd215/IL-15R α Monoclonal Antibody For If/Icc, Fc - Q60819

Item Number: CM0005805



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD215/IL-15R α , validated for IF/ICC and FC. Detects the high-affinity IL-15 receptor, crucial in immune cell signaling.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD215/IL-15R α mAb
Remarks/Alias	IL15ra
Species	Mouse
GeneID (Human)	3601
GeneID	16169
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-205 of mouse CD215/IL-15R α (NP_032384.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Mouse
SWISS	Q60819
Protein Weight	28kDa
Shipping	Ice bag