


Anti-Cd22 Polyclonal Antibody For Wb And Elisa - P20273

Item Number: CM0020703



Introduction

CD22 rabbit polyclonal antibody validated for WB and ELISA. Detects human CD22 (95 kDa), an inhibitory B-cell receptor regulating B-cell signaling and apoptosis. - P20273

[Learn More](#)

Parameter	Value
Product Name	CD22 Rabbit pAb
Remarks/Alias	SIGLEC2, SIGLEC-2, CD22
Species	Human
GeneID (Human)	933
GeneID	933
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CD22 (NP_001762.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P20273
Protein Weight	95kDa
Shipping	Ice bag