


Anti-Cd22 Polyclonal Antibody For Wb, If/Icc, Elisa - P20273

Item Number: CM0030192



Introduction

CD22 Rabbit Polyclonal Antibody, validated for WB, IF/ICC, ELISA. Cross-reacts with Human, Mouse, Rat. Targets B-cell receptor CD22 (SIGLEC2) for B-cell signaling and adhesion research.

[Learn More](#)

Parameter	Value
Product Name	CD22 Rabbit pAb
Remarks/Alias	SIGLEC2; SIGLEC-2; CD22
Species	Human
Gene ID (Human)	933
Gene ID	933
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 708-847 of human CD22 (NP_001762.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20273
Protein Weight	95kDa
Shipping	Ice bag