


Anti-Cd22 Rabbit Polymab® For Wb, If/Icc, If-P, Ihc-P, Elisa - P20273 / P35329

Item Number: CM0007186



Introduction

Recombinant rabbit monoclonal antibody targeting human/mouse CD22 (SIGLEC2), validated for WB, IF, IHC, and ELISA. Ideal for B-cell biology and BCR signaling studies.

[Learn More](#)

Parameter	Value
Product Name	CD22 Rabbit PolymAb®
Remarks/Alias	SIGLEC2; SIGLEC-2
Species	Human/Mouse
GeneID (human)	933
GeneID	933/12483
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P20273/P35329
Protein Weight	95kDa/97kDa
Shipping	Ice bag