

Anti-Cd225/Ifitm1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P13164

Item Number: CM0029311

Introduction

Rabbit polyclonal antibody targeting human CD225/IFITM1 (interferon-induced transmembrane protein 1) for WB, IHC-P, ELISA. Cross-reacts with human, mouse. Ideal for antiviral immunity and cancer studies.

[Learn More](#)

Parameter	Specification
Product Name	CD225/IFITM1 Rabbit pAb
Remarks/Alias	9-27; CD225; IFI17; LEU13; DSPA2a; CD225/IFITM1
Species	Human
GeneID (Human)	8519
GeneID	8519
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-125 of human CD225/IFITM1 (NP_003632.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P13164
Protein Weight	14kDa
Shipping	Ice bag