

Anti-Cd225/Ifitm1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P13164

Item Number: CM0003928

Introduction

High-quality rabbit monoclonal antibody specific to human CD225/IFITM1, validated for WB, IF-P, IHC-P, and ELISA. Detects a 14 kDa interferon-induced transmembrane protein involved in antiviral defense and tumorigenesis. Ideal for innate immunity and oncology studies.

[Learn More](#)

Parameter	Value
Product Name	CD225/IFITM1 Rabbit mAb
Remarks/Alias	9-27; CD225; IFI17; LEU13; DSPA2a; CD225/IFITM1
Species	Human
Gene ID (Human)	8519
Gene ID	8519
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human CD225/IFITM1 (NP_003632.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P13164
Protein Weight	14kDa
Shipping	Ice bag