


Anti-Cd226/Dnam-1 Rabbit Monoclonal Antibody For Wb, Elisa - Q15762

Item Number: CM0005183



Introduction

High-quality anti-CD226/DNAM-1 rabbit monoclonal antibody validated for Western blot (WB) and ELISA applications. Specific for human CD226, a key immune receptor involved in T-cell and NK cell cytotoxicity. Ideal for immunology research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	CD226/DNAM-1 Rabbit mAb
Remarks/Alias	PTA1; DNAM1; DNAM-1; TLISA1; CD226/DNAM-1
Species	Human
Gene ID (Human)	10666
Gene ID	10666
Immunogen	Recombinant protein of human CD226/DNAM-1.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q15762
Protein Weight	39kDa
Shipping	Ice bag