


Anti-Cd23 Polyclonal Antibody For Wb, Ihc-P, Elisa - P06734

Item Number: CM0019619



Introduction

Validate this rabbit polyclonal anti-human CD23 antibody in WB, IHC-P, and ELISA. Detects both human and mouse low-affinity IgE receptor, key in allergy and B cell research.

[Learn More](#)

Parameter	Value
Product Name	CD23 Rabbit pAb
Remarks/Alias	CD23; FCE2; CD23A; IGEBF; CLEC4J; FCER1I; BLAST-2
Species	Human
GeneID (Human)	2208
GeneID	2208
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 48-321 of human CD23 (NP_001993.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P06734
Protein Weight	36kDa
Shipping	Ice bag