

Anti-Cd23 Rabbit Monoclonal Antibody For Wb And Elisa - P06734

Item Number: CM0005747



Introduction

Recombinant rabbit monoclonal antibody against human CD23 (FCER2), validated for WB and ELISA. Targets low-affinity IgE receptor; ideal for B-cell differentiation, allergy, and parasite immunity research.

[Learn More](#)

Parameter	Value
Product Name	CD23 Rabbit mAb
Remarks/Alias	CD23; FCE2; CD23A; IGEBF; CLEC4J; FCERII; BLAST-2
Species	Human
GeneID (human)	2208
GeneID	2208
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 48-321 of human CD23 (NP_001993.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P06734
Protein Weight	36kDa
Shipping	Ice bag