

Anti-Cd239/Bcam Monoclonal Antibody For Wb, Fc, Elisa - P50895

Item Number: CM0006835



Introduction

High-quality rabbit monoclonal antibody against human CD239/BCAM (Basal cell adhesion molecule). Validated for WB, FC, and ELISA applications. Sourced from CamelBio for diagnostic and research use.

[Learn More](#)

Parameter	Value
Product Name	CD239/BCAM Rabbit mAb
Remarks/Alias	AU; LU; CD239; MSK19; CD239/BCAM
Species	Human
GeneID (Human)	4059
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-547 of human CD239/BCAM (NP_005572.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	P50895
Protein Weight	67kDa
Shipping	Ice bag