

## Anti-Cd239/Bcam Polyclonal Antibody For Wb, Elisa - P50895

Item Number: CM0014268



### Introduction

Rabbit polyclonal antibody targeting CD239/BCAM (SWISS: P50895) for WB and ELISA applications. Validated in human, cross-reacts with mouse. Suitable for cell adhesion, migration, and invasion studies.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CD239/BCAM Rabbit pAb                                                                                                   |
| Remarks/Alias    | AU; LU; CD239; MSK19; CD239/BCAM                                                                                        |
| Species          | Human                                                                                                                   |
| Gene ID (Human)  | 4059                                                                                                                    |
| Gene ID          | 4059                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 352-547 of human CD239/BCAM (NP_005572.2) |
| Source           | Rabbit                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                   |
| Application      | WB, ELISA                                                                                                               |
| Cross Reactivity | Mouse                                                                                                                   |
| SWISS            | P50895                                                                                                                  |
| Protein Weight   | 67kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                 |