

Anti-Cd239/Bcam Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P50895

Item Number: CM0006794



Introduction

High-quality rabbit monoclonal antibody targeting human CD239/BCAM. Validated for WB, IHC-P, ELISA. Suitable for cell adhesion and laminin receptor studies. Cross-reacts with human. Protein weight: 67kDa.

[Learn More](#)

Parameter	Value
Product Name	CD239/BCAM Rabbit mAb
Remarks/Alias	BCAM; AU; CD239; LU; MSK19; basal cell adhesion molecule; CD239/BCAM
Species	Human
GeneID	4059
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P50895
Protein Weight	67kDa
Shipping	Ice bag