

Anti-Cd239/Bcam Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P50895

Item Number: CM0002964



Introduction

High-quality rabbit monoclonal antibody targeting human CD239/BCAM (basal cell adhesion molecule). Validated for WB, IHC-P, ELISA with cross-reactivity to mouse. Suitable for cell adhesion, migration, and blood group antigen research.

[Learn More](#)

Parameter	Value
Product Name	CD239/BCAM Rabbit mAb
Remarks/Alias	AU; LU; CD239; MSK19; CD239/BCAM
Species	Human
Gene ID	4059
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 529-628 of human CD239/BCAM (P50895).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P50895
Protein Weight	67kDa
Shipping	Ice bag