

Anti-Cd24 Polyclonal Antibody For Wb, Ihc-P, Elisa - P25063

Item Number: CM0020696

Introduction

CD24 rabbit polyclonal antibody validated for WB, IHC-P, and ELISA; detects human, mouse, and rat CD24. Suitable for studying B-cell differentiation, cancer stem cell markers, and immune modulation research.

[Learn More](#)

Parameter	Value
Product Name	CD24 Rabbit pAb
Remarks/Alias	CD24A; CD24
Species	Human
GeneID (Human)	100133941
GeneID	100133941
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human CD24 (NP_037362.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25063
Protein Weight	8kDa
Shipping	Ice bag