

## Anti-Cd24 Rabbit Monoclonal Antibody For Wb, If, Fc, Elisa - P24807

Item Number: CM0007821

### Introduction

High specificity Anti-CD24 Rabbit Monoclonal Antibody validated for WB, IF/ICC, FC, ELISA in Mouse samples. Recombinant immunogen targeting amino acids 27-53 of mouse CD24. For research use only.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | CD24 Rabbit mAb                                                                                                  |
| Remarks/Alias    | HSA; Cd24; Ly-52; nectadrin                                                                                      |
| Species          | Mouse                                                                                                            |
| GenelD (Human)   | 100133941                                                                                                        |
| GenelD           | 12484                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 27-53 of mouse CD24 (NP_033976.1). |
| Source           | Rabbit                                                                                                           |
| Category         | Monoclonal Antibodies                                                                                            |
| Application      | WB, IF/ICC, IF-P, FC, ELISA                                                                                      |
| Cross Reactivity | Mouse                                                                                                            |
| SWISS            | P24807                                                                                                           |
| Protein Weight   | 8 kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |