

Anti-Cd24 Rabbit Polymab® For Wb, If/Icc, If-P, Elisa - P25063 / P24807

Item Number: CM0008940

Introduction

High-quality anti-CD24 rabbit recombinant monoclonal antibody. Validated for WB, IF/ICC, IF-P, ELISA. Cross-reacts with human and mouse CD24. Useful for studies of B-cell activation, autoimmunity, and cancer.

[Learn More](#)

Parameter	Value
Product Name	CD24 Rabbit PolymAb®
Remarks/Alias	CD24A
Species	Human/Mouse
Gene ID (Human)	100133941
Gene ID (Mouse)	12484
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human/mouse CD24 (NP_037362.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P25063 / P24807
Protein Weight	8kDa / 13kDa
Shipping	Ice bag