

Anti-Cd244 Rabbit Monoclonal Antibody For Wb, Fc, Elisa - Q9Bzw8

Item Number: CM0007598

Introduction

High-quality CD244 rabbit monoclonal antibody validated for WB, FC, and ELISA. Detects human natural killer cell receptor 2B4 (SLAMF4), an immune modulator. Ideal for immunology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CD244 Rabbit mAb
Remarks/Alias	CD244; 2B4; NAIL; NKR2B4; Nmrk; SLAMF4; CD244 molecule
Species	Human
Gene ID (Human)	51744
Gene ID	51744
Immunogen	Recombinant protein corresponding to a sequence within amino acids 22-221 of human CD244 (NP_001160135.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	Q9BZW8
Protein Weight	30kDa/36kDa/41kDa
Shipping	Ice bag