

Anti-Cd244 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bzw8

Item Number: CM0022245

Introduction

CD244 Rabbit Polyclonal Antibody validated for WB, ELISA with cross-reactivity to human and mouse. Targets natural killer cell receptor 2B4 (SLAMF4); immunogen covers amino acids 22-224. For immunological research.

[Learn More](#)

Parameter	Value
Product Name	CD244 Rabbit pAb
Remarks/Alias	2B4; NAIL; Nmrk; NKR2B4; SLAMF4; CD244
Species	Human
Gene ID (Human)	51744
Gene ID	51744
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-224 of human CD244 (NP_057466.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BZW8
Protein Weight	42kDa
Shipping	Ice bag