

Anti-Cd248 Polyclonal Antibody For Wb And Elisa - Q9Hcu0

Item Number: CM0017853

Introduction

Affinity-purified rabbit polyclonal antibody against human CD248 (endosialin, TEM1) for Western blot and ELISA. Validated in human, mouse, rat. Research use for angiogenesis, immune response, and tumor biology studies.

[Learn More](#)

Parameter	Value
Product Name	CD248 Rabbit pAb
Remarks/Alias	TEM1; CD164L1; CD248
Species	Human
GeneID (Human)	57124
GeneID	57124
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 18-150 of human CD248 (NP_065137.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HCU0
Protein Weight	81kDa
Shipping	Ice bag