

Anti-Cd25 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P01589

Item Number: CM0005066



Introduction

High-quality rabbit monoclonal antibody targeting human CD25 (IL2RA), validated for WB, IF/ICC, and ELISA. Ideal for studying immune tolerance, regulatory T-cell biology, and interleukin-2 signaling. Recombinant immunogen, MW 31 kDa.

[Learn More](#)

Parameter	Value
Product Name	CD25 Rabbit mAb
Remarks/Alias	p55; CD25; IL2R; IMD41; TCGFR; IDDM10
Species	Human
GeneID	3559
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-213 of human CD25 (NP_000408.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P01589
Protein Weight	31kDa
Shipping	Ice bag