

Anti-Cd25 Rabbit Polyclonal Antibody - P01589

Item Number: CM0028399



Introduction

High-quality rabbit polyclonal antibody against human CD25 (IL-2 receptor alpha). Validated for WB, IF-P, IHC-P, FC, ELISA with cross-reactivity to mouse and rat. Ideal for immune tolerance and Treg research.

[Learn More](#)

Parameter	Value
Product Name	CD25 Rabbit pAb
Remarks/Alias	p55; CD25; IL2R; IMD41; TCGFR; IDDM10
Species	Human
GeneID (Human)	3559
GeneID	3559
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 101-200 of human CD25 (NP_000408.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, FC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01589
Protein Weight	31 kDa
Shipping	Ice bag