


Anti-Cd252 (Ox40 Ligand) Rabbit Monoclonal Antibody For If/Icc, Fc, Elisa - P43488

Item Number: CM0012493



Introduction

Rabbit monoclonal antibody against mouse CD252 (OX40 Ligand) targeting amino acids 51-198. Validated for IF/ICC, FC, and ELISA. Suitable for T-cell co-stimulation studies. Cross-reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	CD252 (OX40 Ligand) Rabbit mAb
Remarks/Alias	Ath1; gp34; Ath-1; Ox40l; TXGP1; CD134L; OX-40L; Tnlg2b; Txgp1l
Species	Mouse
GeneID (Human)	7292
GeneID	22164
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC, ELISA
Cross Reactivity	Mouse
SWISS	P43488
Protein Weight	22kDa
Shipping	Ice bag