


Anti-Cd262/Dr5/Trailr2 Monoclonal Antibody For Fc, Elisa - O14763

Item Number: CM0004962



Introduction

High-performance rabbit monoclonal antibody specific for human death receptor 5 (DR5, TRAILR2, CD262). Optimized for flow cytometry and ELISA, supporting apoptosis and oncology research. O14763.

[Learn More](#)

Parameter	Value
Product Name	CD262/DR5/TRAILR2 Rabbit mAb
Remarks/Alias	TNFRSF10B; CD262; DR5; KILLER; KILLER/DR5; TRAIL-R2; TRAILR2; TRICK2; TRICK2A; TRICK2B; TRICKB; ZTNFR9; TNF receptor superfamily member 10b; CD262/DR5/TRAILR2
Species	Human
GeneID (human)	8795
GeneID	8795
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 56-182 of human CD262/DR5/TRAILR2 (NP_003833.4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, ELISA
Cross Reactivity	Human
SWISS	O14763
Protein Weight	12kDa/45kDa/47kDa
Shipping	Ice bag