

Anti-Cd26/Dpp4 (Dppiv) Monoclonal Antibody For Flow Cytometry - P27487

Item Number: CM0007734



Introduction

Recombinant Rabbit monoclonal antibody against human CD26/DPP4 (DPPIV, ADCP2). Validated for flow cytometry. Suitable for T-cell activation studies and MERS-CoV receptor research. UniProt P27487.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD26/DPP4 mAb
Remarks/Alias	CD26; ADABP; ADCP2; DPPIV; TP103
Species	Human
GeneID (Human)	1803
GeneID	1803
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 29-766 of human CD26/DPP4 (NP_001926.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P27487
Protein Weight	88kDa
Shipping	Ice bag