


Anti-Cd26/Dpp4 Polyclonal Antibody For Wb, Elisa - P27487

Item Number: CM0020984



Introduction

High-quality rabbit polyclonal antibody against human CD26/DPP4, validated for WB and ELISA with rat cross-reactivity. Ideal for T-cell activation and dipeptidyl peptidase research.

[Learn More](#)

Parameter	Value
Product Name	CD26/DPP4 Rabbit pAb
Remarks/Alias	CD26; ADABP; ADCP2; DPPIV; TP103; CD26/DPP4
Species	Human
GeneID (Human)	1803
GeneID	1803
Immunogen	Recombinant protein corresponding to a sequence within amino acids 29-766 of human CD26/DPP4 (NP_001926.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P27487
Protein Weight	88kDa
Shipping	Ice bag