

Anti-Cd26/Dpp4 Rabbit Monoclonal Antibody For Wb And Elisa - P27487

Item Number: CM0010386



Introduction

CD26/DPP4 rabbit monoclonal antibody detects endogenous human DPP4, a T-cell activation costimulatory receptor and serine protease, applicable in WB and ELISA with cross-reactivity to mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	CD26/DPP4 Rabbit mAb
Remarks/Alias	CD26; ADABP; ADCP2; DPPIV; TP103; CD26/DPP4
Species	Human
Gene ID (Human)	1803
Gene ID	1803
Immunogen	Synthetic peptide corresponding to amino acids 667-766 of human CD26/DPP4 (P27487)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P27487
Protein Weight	88kDa
Shipping	Ice bag