

Anti-Cd26/Dpp4 Rabbit Monoclonal Antibody For Wb, Ip, Fc, Elisa - P27487

Item Number: CM0006933



Introduction

High-quality rabbit monoclonal antibody against human CD26/DPP4 (DPP4), validated in WB, IP, FC, ELISA. Recognizes the 88kDa protein and is suitable for T-cell activation, protease, and viral receptor studies.

[Learn More](#)

Parameter	Specification
Product Name	CD26/DPP4 Rabbit mAb
Remarks/Alias	CD26; ADABP; ADCP2; DPPIV; TP103; CD26/DPP4
Species	Human
Gene ID	1803
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-766 of human CD26/DPP4 (NP_001926.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, FC, ELISA
Cross Reactivity	Human
SWISS	P27487
Protein Weight	88kDa
Shipping	Ice bag