

Anti-Cd26/Dpp4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P27487

Item Number: CM0013006



Introduction

High-quality rabbit polyclonal antibody against human CD26/DPP4. Validated for WB, IHC-P, ELISA, cross-reacts with human, mouse, rat. Ideal for T-cell activation and dipeptidyl peptidase studies.

[Learn More](#)

Parameter	Value
Product Name	CD26/DPP4 Rabbit pAb
Remarks/Alias	CD26; ADABP; ADCP2; DPPIV; TP103; CD26/DPP4
Species	Human
GeneID (Human)	1803
GeneID	1803
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-766 of human CD26/DPP4 (NP_001926.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27487
Protein Weight	88kDa
Shipping	Ice bag