

Anti-Cd27 Monoclonal Antibody For Wb, If-P, Elisa - P26842

Item Number: CM0003832



Introduction

CD27 Rabbit Monoclonal Antibody validated for Western blot (WB), immunofluorescence (IF-P), and ELISA. Recognizes human and rat CD27, a 29 kDa TNF receptor involved in T-cell costimulation. Suitable for immunology and oncology research.

[Learn More](#)

Parameter	Specification
Product Name	CD27 Rabbit mAb
Remarks/Alias	T14; S152; Tp55; TNFRSF7; S152. LPFS2; CD27
Species	Human
Gene ID (Human)	939
Gene ID	939
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CD27 (P26842).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P26842
Protein Weight	29kDa
Shipping	Ice bag