

Anti-Cd27 Polyclonal Antibody For Wb, Elisa - P26842

Item Number: CM0018396



Introduction

High-quality rabbit anti-CD27 polyclonal antibody validated for Western blot and ELISA, detecting human CD27 (~29 kDa) with mouse cross-reactivity. Ideal for immune checkpoint and T-cell research.

[Learn More](#)

Parameter	Value
Product Name	CD27 Rabbit pAb
Remarks/Alias	T14; S152; Tp55; TNFRSF7; S152. LPFS2; CD27
Species	Human
GeneID (Human)	939
GeneID	939
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P26842
Protein Weight	29kDa
Shipping	Ice bag