

Anti-Cd27 Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P26842

Item Number: CM0022505



Introduction

Rabbit polyclonal CD27 antibody validated for WB, IF/ICC, IHC-P, ELISA in human, mouse, rat. Targets CD27 antigen (TNFRSF7), a costimulatory immune-checkpoint receptor involved in T-cell activation and memory. Recombinant protein immunogen.

[Learn More](#)

Parameter	Value
Product Name	CD27 Rabbit pAb
Remarks/Alias	T14; S152; Tp55; TNFRSF7; S152. LPFS2; CD27
Species	Human
Gene ID (Human)	939
Gene ID	939
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-191 of human CD27 (NP_001233.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26842
Protein Weight	29kDa
Shipping	Ice bag