

Anti-Cd271 (P75Ntr) Rabbit Polyclonal Antibody For Ihc-P And Elisa - Q9Z0W1

Item Number: CM0021374

Introduction

Rabbit polyclonal antibody against mouse CD271 (p75NTR, TNFRSF16) validated for IHC-P and ELISA. Cross-reacts with human. Suitable for neurobiology, apoptosis, and rabies virus entry research.

[Learn More](#)

Parameter	Value
Product Name	CD271 Rabbit pAb
Remarks/Alias	p75; LNGFR; p75NTR; p75NGFR; Tnfrsf16
Species	Mouse
GeneID (Human)	4804
GeneID	18053
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-254 of mouse CD271 (NP_150086.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q9Z0W1
Protein Weight	46kDa
Shipping	Ice bag